

HEARING AID MAINTENANCE

Keep Your Hearing Aids in Top Shape!

Daily Care



- Wipe with microfiber cloth
- Brush off earwax (mic & receiver)
- Open battery door overnight

Weekly Check




- Replace wax guards
- Inspect BTE tubing (moisture/cracks)
- Clear microphone ports

Avoid These!



- Moisture (showers, rain, sweat)
- Extreme temps (hot cars, freezing)
- Hairspray/perfume (apply first)

Monthly Deep Clean



- Use cleaning kit (brush/pick)
- Wash BTE earmolds (soap + water)
- Consider professional cleaning

Pro Tips



- Store in drying kit overnight
- Use protective case when traveling
- Keep pets/kids away (they chew!)

Quick Troubleshooting



<u>Problem</u>	<u>Fix</u>
No sound	Replace battery
Whistling	Reinsert/recheck fit
Muffled	Clean wax guards

Did You Know?

Proper care adds 2-3+ years to hearing aids lifespan!

